Rev. DL028J010 / DL028E000

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## Safety Precautions

To help protect you and/or your property from potential damage and ensure personal safety, please read this IFU thoroughly before using the product.

### [Notational Convention]

Notation	Explanation
 CAUTION	Indicates a hazard with a low level of risk which, if not avoided, could result in minor or moderate injury.

#### CAUTION

#### ■ Use only in well ventilated areas

In case of insufficient ventilation, flammable and toxic solvents can cause fire, explosion, or poisoning.

#### ■ First Aid

Skin exposure

Wash exposed area with plenty of soap and water.

Eye exposure

Open eyes as wide as possible and wash with clean water for at least 15 minutes. Immediately call for medical attention.

Ingestion

Please wash mouth with excess water and immediately call for medical attention.

#### ■ Do not spill solvents

Spillage and leakage can cause fire, electric shock, poisoning, injury, and corrosion.

Wear appropriate protective gear when cleaning up a spill.

#### ■ Wear eye protection and protective gloves

Organic solvents and acids are harmful and should not come in direct contact with the skin.

#### ■ Handle package with care

Inappropriate handling may cause rupturing and/or spattering of the product.

#### ■ Only use this product as intended

This product is intended for *in vitro* diagnostic use for the measurement of HbA<sub>1c</sub> in blood specimens.

#### ■ Proper disposal

Dispose in accordance with local laws and regulations.

#### NOTE

Keep this IFU with the product for future reference.

Европейски потребители / Evropské zákazníci / Europæiske kunder / Europäische Kunden / European customers / Ευρωπαϊκοί πελάτες / Clientes europeos / Europejski klienci / Consommateurs européens / Európai ügyfelek / Clienti europei / Eiropas klienti / Klienti europejscy / Clientes europeos / Clienti europei / Evropski kupci / Európski zákazníci / Europeiska kunder

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European Conformity

Manufacturer

Authorized  
representative in the  
European Community

Catalogue number  
/ Part number

*In vitro* diagnostic  
medical device

Consult instructions for use


Use by date

Batch code / Lot number

Temperature limitation

For specified  
column lot only

Serial number /  
Column number

Net volume  
(after reconstitution for  
lyophilized material)

## CONTENTS

1. Introduction .....	18
2. Prior to Use .....	19
3. Warnings and Precautions .....	19
4. Content .....	19
5. Related Components .....	20
6. Storage and Stability .....	20
7. Specimens .....	20
8. Assay Principle .....	20
9. Installation .....	21
10. Assay Procedures .....	22
11. Precautions for Use .....	23
12. Reference Data .....	25
13. Measurement Values .....	25
14. Evaluation of Results .....	25
15. Expected Values .....	26
16. Performance Characteristics .....	27
17. Precision .....	29
18. Interferences .....	29
19. Specimen Collection and Handling .....	31

## 1. Introduction

The GX Assay Kit is designed exclusively for use with the Tosoh Automated Glycohemoglobin Analyzer HLC-723GX (referred to as HLC-723GX in the IFU), which is based on the principle of High-Performance Liquid Chromatography assay. It is not designed for and should never be used with any other type of system.

The GX Assay Kit is intended for *in vitro* diagnostic use for the measurement of Hemoglobin A<sub>1c</sub> (HbA<sub>1c</sub>) in blood specimens. Hemoglobin A<sub>1c</sub> measurements are used in the clinical management of diabetes to assess the long-term efficacy of diabetic control.

Glycohemoglobin (GHb) is a general term for complexes where the whole blood glucose is non-enzymatically bound to the  $\alpha$  or  $\beta$  -chains of human hemoglobin. Within these complexes, HbA<sub>1c</sub>, which is a complex of glucose and the N-terminus of the  $\beta$  -chain, is the most prevalent.

HbA<sub>1c</sub> is non-enzymatically synthesized in two steps.

Step 1: The glucose aldehyde group reacts with the free amino group on the valine at the N-terminus of the  $\beta$  -chain to form the Schiff base (labile HbA<sub>1c</sub> (L-A<sub>1c</sub>)).

Step 2: Stable ketoamine (stable HbA<sub>1c</sub> (s-A<sub>1c</sub>)) is then formed by a reaction known as the amadori rearrangement.

Due to the fact that L-A<sub>1c</sub> (the intermediate product of this reaction) changes rapidly in response to the changes in whole blood glucose concentrations, s-A<sub>1c</sub> is now generally used to measure HbA<sub>1c</sub>. It provides the best indication of average glucose levels over the most recent 1 to 3 month period because it does not fluctuate in response to physiological factors.

In the past it was difficult to achieve quick separation of L-A<sub>1c</sub> and s-A<sub>1c</sub> on a column, therefore sample pretreatment to remove L-A<sub>1c</sub> was necessary before assaying s-A<sub>1c</sub>. The HLC-723GX is featured by quick separation of L-A<sub>1c</sub> and s-A<sub>1c</sub> on the column (2.2 min), with no need for sample treatment when assaying s-A<sub>1c</sub>.

The hemoglobin found in a normal adult is hemoglobin A, which comprises about 97 % of the total hemoglobin. On the other hand, various hemoglobin variants have also been reported in different regions, especially in areas with a large immigrant population.

Abnormal hemoglobins or hemoglobin variants are eluted out together with other Hb components and can lead to misinterpretation of s-A<sub>1c</sub> values. The HLC-723GX with the TSKgel GX and the GX Elution Buffers will not be affected by most commonly-known hemoglobin variants in reporting s-A<sub>1c</sub> values.

## 2. Prior to Use

Inspect the packaging and the exterior of the aluminum pack for any signs of damage prior to use. If any damages are visible, contact your local Tosoh sales representative.

Confirm that the following documents are included in the package.

- Instructions For Use      1 copy

## 3. Warnings and Precautions

- 1) This product is for *in vitro* diagnostic use only.
- 2) This product is intended for use on the Tosoh Automated Glycohemoglobin Analyzer HLC-723GX only.
- 3) Do not use this product beyond the expiration date.

## 4. Content

The GX Assay Kit contains the listed products below.

Catalogue No.	Description	Package contents	
0023161	GX Assay Kit	GX Elution Buffer No. 1	1×800 mL
		GX Elution Buffer No. 2	1×800 mL
		GX Elution Buffer No. 3	1×800 mL
		GX Hemolysis & Wash Solution	2×2000 mL
		Printer Paper	2
		Filter Element	2

The following types of elution buffers are designed for exclusive use with the HLC-723GX.

GX Elution Buffers are organic acid buffers. Each contains less than 0.05 % sodium azide as a preservative.

GX Hemolysis & Wash Solution contains detergent and less than 0.1 % sodium azide as preservative.

Applicable analyzers : Tosoh Automated Glycohemoglobin Analyzer  
HLC-723GX

## 5. Related Components

	Catalogue No.
Hemoglobin A1c Calibrator Set	0018767
Hemoglobin A1c Control	0019405
Hemoglobin A1c Control Set	0021974
TSKgel GX	0023160

## 6. Storage and Stability

GX Elution Buffer (No. 1, No. 2, No. 3)

All unopened materials are stable until the expiration date on the label when stored at 4 to 30 °C. The expiration date is indicated on the package and aluminum pack label.

Stable for 4 months after opening when stored at 4 to 30 °C.

GX Hemolysis & Wash Solution

All unopened materials are stable until the expiration date on the label when stored at 4 to 30 °C. The expiration date is listed on the bottle label.

Stable for 4 months after opening when stored at 4 to 30 °C.

## 7. Specimens

Whole blood samples in primary tubes containing EDTA.

## 8. Assay Principle

The HLC-723GX is based on the High-Performance Liquid Chromatography (HPLC) principle with a cation exchange column employing the differences in ionic interactions between hemoglobin components to separate them in a span of 2.2 min.

A step gradient elution is used to separate HbF, s-A<sub>1c</sub>, Total A<sub>1</sub>, and Hb variants with a set of Elution Buffers (GX Elution Buffer No. 1, 2, and 3) of different salt concentrations.

## 9. Installation

Be sure to read the IFU included with the TSKgel GX, GX Assay Kit, Hemoglobin A1c Calibrator Set, Hemoglobin A1c Control, Hemoglobin A1c Control Set as well as HLC-723 GX Operator's Manual.

### Installation of GX Elution Buffer

- 1) Press the STOP key on the HLC-723GX system operation panel to enable the STAND-BY state.
- 2) Remove the cap from the elution buffer pack.
- 3) Make sure to insert the tube matching the color of the elution buffer pack being installed.
- 4) Gently squeeze the elution buffer pack by hand to remove all excess air, and tighten cap firmly into place to create a vacuum that will prevent air from entering during operation.
- 5) Press the MAINT key and then the REAGENT CHANGE key.
- 6) Press the keys of the reagents to be replaced on the REAGENT CHANGE screen.
- 7) Press the CHANGE key. A confirmation message will be displayed.

If everything is all right, press the OK key.

### Installation of GX Hemolysis & Wash Solution

- 1) Press the STOP key on the HLC-723GX system operation panel to enable the STAND-BY state.
- 2) Remove the cap from the GX Hemolysis & Wash Solution bottle.
- 3) Insert the tube into the new bottle and screw cap back on.
- 4) Press the MAINT key and then the REAGENT CHANGE key.
- 5) Press the H/W key on the REAGENT CHANGE screen.
- 6) Press the CHANGE key. A confirmation message will be displayed.

If everything is all right, press the OK key.

### Installation of Printer Paper

- 1) Lift the printer cover (upper lid) to open and push the paper holding lever down to the very front and wrap the remaining paper onto the roll.
- 2) Lift the roll up and remove the mandrel.
- 3) Insert the mandrel into the new roll and set it with great care for the direction.
- 4) Return the paper holding lever to the very back and insert the paper into the printer. The paper will automatically be fed.

### Installation of Filter Element

- 1) Confirm that analyzer system eluent delivery has been terminated (STAND-BY is displayed on the Main screen (first screen)). Then, remove the filter outlet (top) flow line.
- 2) Turn the top of the filter holder assembly by turning it counterclockwise and remove the filter holder by pulling it straight out.
- 3) Lightly press the top of the holder to take out the used filter element.
- 4) Place the new element with great care of surface of filter. The gray colored surface is the outlet (up) side.
- 5) Firmly tighten the top of the filter holder assembly by hand until no further tightening is possible.
- 6) Press the arrow key at the bottom right of the Main screen (first screen) to show the Main screen (second screen). Press the pump key on the Main screen (second screen) to start the pump and check the eluent flows through the filter outlet (top). Then, stop the pump by pressing the pump key again. Take special care that eluent coming from the tube does not spill onto the analyzer unit by preparing wiping paper beforehand.
- 7) Connect the outlet side flow tube.

## 10. Assay Procedures


Refer to the HLC-723GX Operator's Manual for detailed instructions.

## 11. Precautions for Use

- 1) Carefully read the instructions contained in this IFU and related Instructions For Use provided with HLC-723GX, TSKgel GX, Hemoglobin A1c Calibrator Set, Hemoglobin A1c Control Set, Hemoglobin A1c Control Set.
- 2) GX Assay Kit is designed exclusively for use in combination with the HLC-723GX and TSKgel GX, never in any other combination.
- 3) For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state and federal regulations.

### **GX Elution Buffer**

- 1) Return the elution buffer to 15 to 30 °C before use.
- 2) Always use the GX Elution Buffer in combination with the TSKgel GX of the identical lot number. The column lot number is indicated by a single uppercase alphabetical character (A, B, etc.) on the label of column box. The Elution Buffer label displays an alphabetic character corresponding to the column lot number, as shown below.

- 3) Always use GX Elution Buffers from the same GX Assay kit number. The GX Assay kit number is indicated by 2 digits and two uppercase alphabetical characters (A01K, A02M, etc.) on the aluminum pack label, as shown below.

**KIT A01K**

- 4) When first opening the aluminum pack to insert the suction tube, gently squeeze the elution buffer pack by hand to remove all excess air, and tighten the cap firmly to create a vacuum that will prevent air from entering during operation.

- 5) Never use reagents that have exceeded the expiration date printed on their labels. Assay results for expired reagents will not be reliable. Also note that reagents must be used within four months after opening (provided that the container is correctly maintained in vacuum state).
- 6) When there is leftover elution buffer in the aluminum pack that must be removed from the analyzer and stored, again gently squeeze the elution buffer pack by hand to remove all excess air, and tighten the cap firmly to create a vacuum that will prevent air from entering during storage and store between 4 and 30 °C. This product is stable for four months after opening.
- 7) Always replace with an aluminum pack when elution buffer is almost empty. Avoid refilling leftover elution buffer into aluminum pack as this can produce unreliable assay results.

#### **GX Hemolysis & Wash Solution**


- 1) Always return GX Hemolysis & Wash Solution to 15 to 30 °C before use.
- 2) Assay results for expired reagents will not be reliable. Also note that reagents must be used within four months after opening (provided that the container is correctly maintained in vacuum state).
- 3) This product is stable for four months after opening when stored at 4 to 30 °C.
- 4) Always replace with a new container when GX Hemolysis & Wash Solution is about to run out. Avoid refilling leftover GX Hemolysis & Wash Solution containers as this can produce unreliable assay results.

#### **Filter Element**


- 1) Place the new filter element with great care of filter surface. The gray colored surface is the outlet (up) side.
- 2) As once the filter element is tightened it is deformed, it cannot be used again.
- 3) The used filter element has been in contact with blood sample. Wear protective clothing (glasses, gloves, mask, etc.) and take sufficient care to prevent potential infection during installation and handling.

## 12. Reference Data (from the HLC-723GX Assay Examples)

Control

Sample

## 13. Measurement Values

Measurement values (%) on the above report, the second column, indicate the percentage of each peak in relation to the Total Area (excluding the front peak (FP)) after calibration. Note the minimum unit of measurement displayed is 0.1 %.

## 14. Evaluation of Results

Quality Control

In order to monitor and evaluate the precision and analytical performance, it is recommended that control samples should be assayed daily.

If one or more control sample value(s) is out of the acceptable

range, it is necessary to investigate the validity of the calibration curve before reporting patient results. Standard laboratory procedures should be followed in accordance to the strictest regulatory agency under which the laboratory operates.

## 15. Expected Values

Reference ranges (nondiabetic):HbA<sub>1c</sub> (s-A<sub>1c</sub>) 4.0 %-6.0 % (mean 5.0 %, SD 0.5 %)

Ref: American Diabetes Association. Standards of medical care for patients with diabetes mellitus. Diabetes Care 2000; 23 (Suppl. 1), S32-42.

The values referred to within this IFU have been determined with a NGSP certified method. It is known that the relationship between HbA<sub>1c</sub> results from the NGSP network (%) and the IFCC network (mmol/mol) is expressed by using the following equation  
(See [http://www.ifcchba1c.net/IFCC\\_08.asp](http://www.ifcchba1c.net/IFCC_08.asp)):

$$\text{NGSP} = 0.09148 \times \text{IFCC} + 2.152$$

## 16. Performance Characteristics

### 1) Dilution Linearity

The tests were performed with Hemoglobin A1c Control diluted with various percentage factors to determine the effect on the s-A<sub>1c</sub> (%) value. Linearity was observed for the Total Area in the 700 to 4500 range, and measurements in the range of 800 to 4000 are recommended to give more reliable results.

Total Area	s-A <sub>1c</sub> (%)
493.53	10.3
986.09	10.4
1492.07	10.5
1933.15	10.5
2427.67	10.5
2907.35	10.5
3403.81	10.6
3803.19	10.6
4292.48	10.6
4698.64	10.6

### 2) Recovery

The table below shows the recovery ratios for two types of samples mixed in different percentage factors.


High s-A <sub>1c</sub> sample (ratio)	Low s-A <sub>1c</sub> sample (ratio)	Measurement value s-A <sub>1c</sub> (%)	Theoretical value s-A <sub>1c</sub> (%)	Recovery (%)
0	10	2.7	-	-
1	9	4.5	4.4	102
2	8	6.2	6.1	102
3	7	7.9	7.8	101
4	6	9.4	9.5	99
5	5	11.3	11.2	101
6	4	13.1	12.9	102
7	3	14.7	14.6	101
8	2	16.3	16.3	100
9	1	18.1	18.0	101
10	0	19.7	-	-


3) Correlation

Correlations (1) between the HLC-723GX and the HLC-723G8 and (2) between the HLC-723GX and the HLC-723G7 are shown in the graphs below.

(1) Correlation between the HLC-723GX and the HLC-723G8.

(2) Correlation between the HLC-723GX and the HLC-723G7.

## 17. Precision

### 1) Intra-assay (n=10)

The intra-assay precision (n=10) of the s-A<sub>1c</sub> (%) for three types of control samples are shown below.

	Low	Medium	High
Mean	5.40	7.99	10.50
SD	0.00	0.03	0.00
CV (%)	0.0	0.4	0.0

### 2) Inter-assay (n=10)

The inter-assay precision (n=10), measured once a day for a 10-day period, of s-A<sub>1c</sub> (%) for two types of control samples are shown below.


	Low	High
Mean	5.34	10.42
SD	0.05	0.06
CV (%)	0.9	0.6

## 18. Interferences


A substance is considered to interfere when the recovery of a known specimen falls outside a 10 % acceptable range.

- 1) No effect was observed up to 10 g/L for glucose, 0.25 g/L for sodium cyanate and 0.25 g/L for acetaldehyde.


Glucose

### Sodium cyanate

### Acetaldehyde

2) No effect was observed for the following agents up to the concentrations indicated in the table.

Free bilirubin	to 0.19 g/L
Conjugated bilirubin	to 0.20 g/L
Chyle	to 1500 Formazine turbidity
Aspirin	to 0.50 g/L
Triglyceride	to 20 g/L

3) The presence of HbD, HbS, and HbC does not interfere with the assay. High HbF may interfere with the assay.

4) HbE may interfere with assay. In the presence of HbE, a shoulder or small peak may be observed behind s-A<sub>1c</sub> peak.

- 5) Since a high number of Hb variants is observed worldwide, care should be taken when interpreting the results. Designation of a particular name cannot exclude the presence of other variants.
- 6) Low concentrations of hemoglobin may cause the Total Area of the chromatogram to be < 700.
- 7) A reduced concentration of HbA<sub>1c</sub> may be observed in patients due to limited life span of red blood cells; interpret clinical results with caution.
- 8) Hemolysis does not interfere with HbA<sub>1c</sub> assay results.

## **19. Specimen Collection and Handling**

A whole blood sample is required for the assay.

No special preparation is necessary. A venous blood sample is collected aseptically. The blood samples in primary tubes containing EDTA, may be stored at 25 °C for 24 hours prior to analysis or at 4 °C for 14 days.

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This manual may not be reprinted or copied in whole or in part without written consent of TOSOH CORPORATION. The contents of the manual are subject to change without notice.

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HLC and TSKgel are registered trademarks of TOSOH CORPORATION.